



Product Name: Recombinant Human MET Protein
Catalog #: PEH3388



Summary

Name	MET_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P08581
Host	E.coli
Species	Human
Molecular Mass	58 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AKSEMNVNMK YQLPNFTAET PIQNVILHEH HIFLGATNYI YVLNEEDLQK VAEYKTGPVL EHPDCFPCQD CSSKANLSGG VWKDNINMAL VVDITYDDQL ISCGSVNRGT CQRHVFPNH TADIQSEVHC IFSPQIEEPS QCPDCVVSAL GAKVLSSVKD RFINFFVGNT INSSYFPDHP LHSISVRRLK ETKDGFMTLT DQSYIDVLPE FRDSYPIKYV HAFESNNFIY FLTQVRETLD AQTFHTRIIR FCSINSLHS YMEMPLECIL TEKRRKRSTK KEVFNILQAA YVSKPGAQLA RQIGASLND ILFGVFAQSK PDSAEPMDRS AMCAFIKYV NDFFNKIVNK NNVRCLQHFY GPNHEHCFNR TLLRNSSGCE ARRDEYRTEF

image



Background

Alternative Names	Hepatocyte growth factor receptor (HGF receptor) (EC 2.7.10.1) (HGF/SF
--------------------------	--

Product Name: Recombinant Human MET Protein
Catalog #: PEH3388



receptor) (Proto-oncogene c-Met) (Scatter factor receptor) (SF receptor)
(Tyrosine-protein kinase Met)

Note

For Research Use Only , Not for Diagnostic Use