


**Product Name: Recombinant Human GSTP1 Protein**  
**Catalog #: PEH3389**



## Summary

|                              |                                                                                                                                                                                                                                                    |
|------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | GSTP1_HUMAN                                                                                                                                                                                                                                        |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                                                                |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                                                                   |
| <b>Tag</b>                   | C-His                                                                                                                                                                                                                                              |
| <b>Accession #</b>           | P09211                                                                                                                                                                                                                                             |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                                                             |
| <b>Species</b>               | Human                                                                                                                                                                                                                                              |
| <b>Molecular Mass</b>        | 24 kDa                                                                                                                                                                                                                                             |
| <b>Buffer</b>                | Phosphate buffered saline                                                                                                                                                                                                                          |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                                                                 |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                                     |
| <b>Sequence</b>              | MPPYTVVYFP VRGRCAALRM LLADQGQSWK EEVVTVETWQ EGSLKASCLY<br>GQLPKFQDGD LTLYQSNTIL RHLGRTLGLY GKDQQEAALV DMVNDGVEDL<br>RCKYISLIYT NYEAGKDDYV KALPGQLKPF ETLLSQNQGG KTFIVGDQIS<br>FADYNLLDLL LIHEVLAPGC LDAFPLLSAY VGRLSARPKL KAFLASPEYV<br>NLPINGNGKQ |

## image



## Background

**Alternative Names** Glutathione S-transferase P (EC 2.5.1.18) (GST class-pi) (GSTP1-1)

## Note

**Product Name: Recombinant Human GSTP1 Protein**  
**Catalog #: PEH3389**



---

For Research Use Only , Not for Diagnostic Use