


Product Name: Recombinant Human CD36 Protein
Catalog #: PEH3391



Summary

Name	CD36_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P16671
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QVVLEEGTIA FKNWVKTGTE VYRQFWIFDV QNPQEVMMNS SNIQVKQRGP YTYRVRFLAK ENVTQDAEDN TVSFLQPNGA IFEPSLSVGT EADNFTVLNL AVAAASHIQ NQFVQMILNS LINKSKSSMF QVRTLRELLW GYRDPFLSLV PYPVTTTVGL FYPYNTADG VYKVFNGKDN ISKVAIIDTY KGKRNLSEWE

image



Background

Alternative Names	Platelet glycoprotein 4 (Fatty acid translocase) (FAT) (Glycoprotein IIIb) (GPIIIB) (Leukocyte differentiation antigen CD36) (PAS IV) (PAS-4) (Platelet collagen receptor) (Platelet glycoprotein IV) (GPIV) (Thrombospondin receptor) (CD antigen CD36)
--------------------------	--

Product Name: Recombinant Human CD36 Protein
Catalog #: PEH3391



Note

For Research Use Only , Not for Diagnostic Use