


Product Name: Recombinant Human ENG Protein
Catalog #: PEH3393



Summary

Name	EGLN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P17813
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TLQASKQNGT WPREVLVLS VNSSVFLHLQ ALGIPLHLAY NSSLVTFQEP PGVNTTELPS FPKTQILEWA AERGPITSAA ELNDPQSILL RLGQAQGSLS FCMLEASQDM GRTLEWRPRT PALVRGCHLE GVAGHKEAHI LRVLPGHSAG PRTVTVKVEL SCAPGDLDAV LILQGPPYVS WLIDANHNMQ IWTTGEYSFK IFPEKNIRGF KLPDTPQGLL

image



Background

Alternative Names Endoglin (CD antigen CD105)

Note

Product Name: Recombinant Human ENG Protein
Catalog #: PEH3393



For Research Use Only , Not for Diagnostic Use