


Product Name: Recombinant Human ERLIN1 Protein
Catalog #: PEH3397



Summary

Name	ERLN1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75477
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EKTKLLIAAQ KQKVVEKEAE TERKKAVIEA EKIAQVAKIR FQQKVMKETS EKRISEIEDA AFLAREKAKA DAEYAAHKY ATSNKHKLTP EYLELKKYQA

image



Background

Alternative Names Erlin-1 (Endoplasmic reticulum lipid raft-associated protein 1) (Protein KE04) (Stomatin-prohibitin-flotillin-HflC/K domain-containing protein 1) (SPFH domain-containing protein 1)

Note

For Research Use Only , Not for Diagnostic Use