


Product Name: Recombinant Human CALB1 Protein
Catalog #: PEH3398



Summary

Name	CALB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P05937
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ITASQFFEIW LHFADGSGY LEGKELQNLI QELQQARKKA GLELSPMKT FVDQYGQRDD GKIGIVELAH VLPTEENFLL LFRCQQLKSC EEFMKTWRKY DTDHSGFIET EELKNFLKDL LEKANKTVDD TKLAEYTDLM LKLFDSNNDG KLELTEMARL LPVQENFLK FQGIKMCCKE FNKAFELYDQ DNGYIDENE LDALLKDLCE KNKQDLINN ITTYKKNIMA LSDGGKLYRT

image



Background

Alternative Names Calbindin (Calbindin D28) (D-28K) (Vitamin D-dependent calcium-binding protein, avian-type)

Note

Product Name: Recombinant Human CALB1 Protein
Catalog #: PEH3398



For Research Use Only , Not for Diagnostic Use