


Product Name: Recombinant Human PAM Protein
Catalog #: PEH3399



Summary

Name	AMD_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P19021
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VKSDMVMME HHKETEYKDK IPLLQPKRE EEEVLDQGDF YSLLSKLLGE REDVVHVHKY NPTEKAESES DLVAEIANVV QKKDLGRSDA REGAEHERGN

image



Background

Alternative Names	Peptidyl-glycine alpha-amidating monooxygenase (PAM) [Includes: Peptidylglycine alpha-hydroxylating monooxygenase (PHM) (EC 1.14.17.3); Peptidyl-alpha-hydroxyglycine alpha-amidating lyase (EC 4.3.2.5) (Peptidylamidoglycolate lyase) (PAL)]
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PAM Protein
Catalog #: PEH3399

