


Product Name: Recombinant Human MSR1 Protein
Catalog #: PEH3400



Summary

Name	MSRE_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P21757
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KNCSVSSTNA NDITQSLTGK GNDSEEEEMRF QEVFMEHMSN MEKRIQHILD MEANLMDTEH FQNFSMTTDQ RFNDILLQLS TLFSSVQGHG NAIDEISKSL ISLNTTLLDL QLNIENTLNGK IQENTFKQQE EISKLEERVY NVSAEIMAMK EEQVHLEQEI KGEVKVLNNI TNDLRLKDWE HSQTLRNITL

image



Background

Alternative Names	Macrophage scavenger receptor types I and II (Macrophage acetylated LDL receptor I and II) (Scavenger receptor class A member 1) (CD antigen CD204)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MSR1 Protein
Catalog #: PEH3400

