

Summary

Name	DCTN1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14203
Host	E.coli
Species	Human
Molecular Mass	68 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LETLRKRAE DKAKLKELEK HKIQLEQVQE WSKMQEQQA DLQRRLEKEAR KEAKEALEAK ERYMEEMADT ADAIEMATLD KEMAEERAES LQQEVEALKE RVDELTTDLE ILKAEIEEKG SDGAASSYQL KQLEEQNARL KDALVRMRDL SSSEKQEHVK LQKLMEKKNQ ELEVVRQRE RLQEELSQAE STIDELKEQV DAALGAEEMV EMLTDRNLNL EEKVRELRET VGDLEAMNEM NDELQENARE TELELREQLD MAGARVREAQ KRVEAAQETV ADYQQTIKKY RQLTAHLQDV NRELTNQQEA SVERQQQPPP ETFDFKIKFA ETKAHAKAIE MELRQMEVAQ ANRHMSLLTA FMPDSFLRPG GDHDCVLVLL LMPRLICKAE LIRKQAQEKF ELSENCSERP GLRGAAGEQL SFAAGLVYSL SLLQATLHRY EHALSQCSVD VYKKVGSLYP EMSAHERSLD FLIELLHKDQ LDETVNVEPL TKAIKYYQHL YSIHAEQPE DCTMQLADHI KFTQSALDCM SVEVGRLRAF LQGGQEATDI ALLLRDLETS CSDIRQFCKK

image



