


Product Name: Recombinant Human ANGPT1 Protein
Catalog #: PEH3405



Summary

Name	ANGP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q15389
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NQRRSPENSG RRYNRIQHGX CAYTFILPEH DGNCRESTTD QYNTNALQRD APHVEPDFSS QKLQHLEHVM ENYTQWLQKL ENYIVENMKS EMAQIQQNAV QNHTATMLEI GTSLLSQTAE QTRKLTDVET QVLNQTSRLE IQLLENSLST YKLEKQLLQQ TNEILKIEK NSLLEHKILE MEGKHKEELD TLKEEKENLQ GLVTRQTYII QLEKQLNRA TTNSVLQKQ QLELMDTVHN LVNLCTKEGV

image



Background

Alternative Names Angiopoietin-1 (ANG-1)

Note

Product Name: Recombinant Human ANGPT1 Protein
Catalog #: PEH3405



For Research Use Only , Not for Diagnostic Use