



Product Name: Recombinant Human ANGPT1 Protein
Catalog #: PEH3406



Summary

Name	ANGP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q15389
Host	E.coli
Species	Human
Molecular Mass	54 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QTRKLTDVET QVLNQTSRLE IQLLENSLST YKLEKQLLQQ TNEILKIHEK NSLLEHKILE MEGKHKEELD TLKEEKENLQ GLVTRQTYII QELEKQLNRA TTNNSVLQKQ QLELMDTVHN LVNLCTKEGV LLKGGKREEE KPFRDCADV QAGFNKSGIY TIYINNMPPEP KKVFCNMDVN GGGWTVIQHR EDGSLDFQRG WKEYKMGFGN PSGEYWLGN FIFAITSRQ YMLRIELMDW EGNRAYSQYD RFHIGNEKQN YRLYLKGHTG TAGKQSSLIL HGADFSTKDA DNDNCMCKCA LMLTGGWWFD ACGPSNLNGM FYTAGQNHGK LNGIKWHYFK GPSYSRSTT

image



Background

Alternative Names Angiopoietin-1 (ANG-1)

Product Name: Recombinant Human ANGPT1 Protein
Catalog #: PEH3406



Note

For Research Use Only , Not for Diagnostic Use