



Product Name: Recombinant Human CLIP2 Protein
Catalog #: PEH3414



Summary

Name	CLIP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UDT6
Host	E.coli
Species	Human
Molecular Mass	67 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSVSSVASSV GGRPSRSGLL TETSSRYARK ISGTTALQEA LKEKQQHIEQ LLAERDLERA EVAKATSHIC EVEKEIALLK AQHEQYVAEA EEKLQRARLL VESVRKEKVD LSNQLEEERR KVEDLQFRVE EESITKGDLE TQTQLEHARI GELEQSLLLE KAQAERLLRE LADNRLTTVA EKSRLVQLEE ELTLRRGEIE ELQQCLLHSG PPPPDHPDAA EILRLRERLL SASKEHQRES GVLRDKYEKA LKAYQAEVDK LRAANEKYAQ EVAGLKDKVQ QATSENMGML DNWWSKLDL ASDHQKSLED LKATLNSGPG AQQKEIGELK AVMEGIKMEH QLELGNLQAK HDLETAMHVK EKEALREKLQ EAQEELAGLQ RHWRAQLEVQ ASQHRLELQE AQDQRRDAEL RVHELEKLDV EYRGQAQAE FLKEQISLAE KKMLDYERLQ RAEAQGKQEV ESLREKLLVA ENRLQAVEAL CSSQHTHMIE SNDISEETIR TKETVEGLQD KLNKRDKVET ALTSQTEMLR AQVSALESKC KSGEKKVDAL LKEKRRLEAE LETVSRKTHD

image



Product Name: Recombinant Human CLIP2 Protein
Catalog #: PEH3414



Background

Alternative Names

CAP-Gly domain-containing linker protein 2 (Cytoplasmic linker protein 115) (CLIP-115) (Cytoplasmic linker protein 2) (Williams-Beuren syndrome chromosomal region 3 protein) (Williams-Beuren syndrome chromosomal region 4 protein)

Note

For Research Use Only , Not for Diagnostic Use