



Product Name: Recombinant Human GALNT2 Protein
Catalog #: PEH3422



Summary

Name	GALT2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q10471
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LAGGAGGGAG RKEDWNEIDP IKKKDLHHSN GEEKAQSMET LPPGKVRWPD FNQEAYVGGT MVRSGQDPYA RNKFNQVESD KLRMDRAIPD

image



Background

Alternative Names	Polypeptide N-acetylgalactosaminyltransferase 2 (EC 2.4.1.41) (Polypeptide GalNAc transferase 2) (GalNAc-T2) (pp-GaNTase 2) (Protein-UDP acetylgalactosaminyltransferase 2) (UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase 2) [Cleaved into: Polypeptide N-acetylgalactosaminyltransferase 2 soluble form]
--------------------------	---

Note

Product Name: Recombinant Human GALNT2 Protein
Catalog #: PEH3422



For Research Use Only , Not for Diagnostic Use