


Product Name: Recombinant Human FCRL5 Protein
Catalog #: PEH3424



Summary

Name	FCRL5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96RD9
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PASDPARSPS DSDSQEPTYH NVPaweELQP VYTNANPRGE NVVYSEVRII QEKKKHAVAS DPRHLRNKGS PIYSEVKVA

image



Background

Alternative Names	Fc receptor-like protein 5 (FcR-like protein 5) (FcRL5) (BXMAS1) (Fc receptor homolog 5) (FcRH5) (Immune receptor translocation-associated protein 2) (CD antigen CD307e)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use