


Product Name: Recombinant Human EIF4B protein
Catalog #: PEH3425



Summary

Name	IF4B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-GST
Accession #	P23588
Host	E.coli
Species	Human
Molecular Mass	52 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QKEQEKLQRQ LDEPKLRRP RERHPSWRSE ETQERERSRT GSESSQTGTS TTSSRNARRR ESEKSLENET LNKEEDCHSP TSKPPKPDQP LKVMPPAPPPK ENAWVKRSSN PPARSQSSDT EQQSPTSGGG KVAPAPQSEE GPGRKDENV DGMNAPKGQT GNSSRGPGDG GNRDHWKESD RKGKKDQDS RSAPEPKKPE ENPASKFSSA SKYAALSVDG EDENEGEDYA

image



Background

Alternative Names Eukaryotic translation initiation factor 4B (eIF-4B)

Note

Product Name: Recombinant Human EIF4B protein
Catalog #: PEH3425



For Research Use Only , Not for Diagnostic Use