


Product Name: Recombinant Human PCDH12 Protein
Catalog #: PEH3427



Summary

Name	PCD12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NPG4
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VLQDTVNLLF NHPRQRNASR ENLNLPEPQP ATGQPRSRPL KVAGSPTGRL AGDQGSEEAP QRPPASSATL RRQRHLNGKV

image



Background

Alternative Names Protocadherin-12 (Vascular cadherin-2) (Vascular endothelial cadherin-2) (VE-cad-2) (VE-cadherin-2) [Cleaved into: Protocadherin-12, secreted form]

Note

For Research Use Only , Not for Diagnostic Use