


Summary

Name	GFRA2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O00451
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DFTENPCLRN AIQAFGNGTD VNVSPKGPSF QATQAPRVEK TPSLPDDLSD STSLGTSVIT TCTSVQEQGL

image



Background

Alternative Names	GDNF family receptor alpha-2 (GDNF receptor alpha-2) (GDNFR-alpha-2) (GFR-alpha-2) (GDNF receptor beta) (GDNFR-beta) (Neurturin receptor alpha) (NRTNR-alpha) (NTNR-alpha) (RET ligand 2) (TGF-beta-related neurotrophic factor receptor 2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GFRA2 Protein
Catalog #: PEH3429

