


Summary

Name	CASQ2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14958
Host	E.coli
Species	Human
Molecular Mass	45 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KDRVVSLSEK NFKQVLKKYD LLCLYYHEPV SSDKVTQKQF QLKEIVLELV AQVLEHKAIG FVMVDAKKEA KLAKKLGFE EGSYLKGD RTIEFDGEFA ADVLVEFLLD LIEDPVEIIS SKLEVQAFER IEDYIKLIGF FKSEDSEYYK AFEEAAEHFQ PYIKFFATFD KGVAKKLSLK MNEVDFYEPF MDEPIAIPNK PYTEEELVEF VKEHQRPTLR RLRPEEMFET WEDDLNGIHI VAFAEKSDPD GYEFLEILKQ VARDNTDNPD LSILWIDPDD FLLVAYWEK TFKIDLFRPQ IGVVNVTDAD SVWMEIPDDD DLPTAELED WIEDVLSGKI NTEDDDEDDD DDDNSDEEDN

image



Background

Alternative Names Calsequestrin-2 (Calsequestrin, cardiac muscle isoform)

Product Name: Recombinant Human CASQ2 Protein
Catalog #: PEH3430



Note

For Research Use Only , Not for Diagnostic Use