


Product Name: Recombinant Human ADAM28 Protein
Catalog #: PEH3435



Summary

Name	ADA28_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UKQ2
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SREKQKKDQR PLSTTGTRPH KQKRKPQMVK AVQPQEMSQM KPHVYDLPVE GNEPPASFHK DTNALPPTVF KDNPVSTPKD

image



Background

Alternative Names	Disintegrin and metalloproteinase domain-containing protein 28 (ADAM 28) (EC 3.4.24.-) (Epididymal metalloproteinase-like, disintegrin-like, and cysteine-rich protein II) (eMDC II) (Metalloproteinase-like, disintegrin-like, and cysteine-rich protein L) (MDC-L)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ADAM28 Protein
Catalog #: PEH3435

