


Product Name: Recombinant Human ITIH4 Protein
Catalog #: PEH3441



Summary

Name	ITIH4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q14624
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HEDSDSALQL QDFYQEVANP LLTAVTFEYP SNAVEEVTQN NFRLLFKGSE MVVAGKLQDR GPDVLTATVS GKLPTQNITF QTESSVAEQE

image



Background

Alternative Names	Inter-alpha-trypsin inhibitor heavy chain H4 (ITI heavy chain H4) (ITI-HC4) (Inter-alpha-inhibitor heavy chain 4) (Inter-alpha-trypsin inhibitor family heavy chain-related protein) (IHRP) (Plasma kallikrein sensitive glycoprotein 120) (Gp120) (PK-120) [Cleaved into: 70 kDa inter-alpha-trypsin inhibitor heavy chain H4; 35 kDa inter-alpha-trypsin inhibitor heavy chain H4]
--------------------------	--

Note

Product Name: Recombinant Human ITIH4 Protein
Catalog #: PEH3441



For Research Use Only , Not for Diagnostic Use