


Product Name: Recombinant Human ATP5IF1 Protein
Catalog #: PEH3443



Summary

Name	ATIF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UII2
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GSDQSENVDR GAGSIREAGG AFGKREQAEE ERYFRAQSRE QLAALKKHHE EEIVHHKKEI ERLQKEIERH KQKIKMLKHD D

image



Background

Alternative Names	ATPase inhibitor, mitochondrial (ATP synthase F1 subunit epsilon) (Inhibitor of F(1)F(o)-ATPase) (IF(1)) (IF1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use