


Product Name: Recombinant Human F10 Protein
Catalog #: PEH3445



Summary

Name	FA10_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P00742
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GHLEECMEE TCSYEEAREV FEDSDKTNEF WNKYKDGQDC ETSPCQNQGK CKDGLGEYTC TCLEGFEGKN

image



Background

Alternative Names	Coagulation factor X (EC 3.4.21.6) (Stuart factor) (Stuart-Prower factor) [Cleaved into: Factor X light chain; Factor X heavy chain; Activated factor Xa heavy chain]
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use