


Product Name: Recombinant Human SEMA4D Protein
Catalog #: PEH3448



Summary

Name	SEM4D_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92854
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PRQCLKFRSA LLIGKKPKS DFCDREQSLK ETLVEPGSFS QQNGEHPKPA LDTGYETE QD TITSKVPTDR EDSQRIDDL S ARDKPFDVKC

image



Background

Alternative Names Semaphorin-4D (A8) (BB18) (GR3) (CD antigen CD100)

Note

For Research Use Only , Not for Diagnostic Use