


Product Name: Recombinant Human XYLT2 Protein
Catalog #: PEH3450



Summary

Name	XYLT2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H1B5
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KPRPLDPGEG SKDTDSSAGR RGSTGRRHGR WRGRAESPGV PVAKVVRAVT SRQRASRRVP PAPPPEAPGR QNLSGAAAGE ALVGAAGFPP HGDTGSVEGA

image



Background

Alternative Names Xylosyltransferase 2 (EC 2.4.2.26) (Peptide O-xylosyltransferase 1) (Xylosyltransferase II) (XT-II) (XylT-II)

Note

For Research Use Only , Not for Diagnostic Use