


Product Name: Recombinant Human LRP1B Protein
Catalog #: PEH3451



Summary

Name	LRP1B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NZR2
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FKCTKNQKCI PVNLRNGQD DCGDEEDERD CPENSCSPDY FQCKTTKHCI SKLWVCDEDP DCADASDEAN CDKKTCPHE FQCKNNNCIP DHWRCDSQND CSDNSDEENC

image



Background

Alternative Names Low-density lipoprotein receptor-related protein 1B (LRP-1B) (Low-density lipoprotein receptor-related protein-deleted in tumor) (LRP-DIT)

Note

For Research Use Only , Not for Diagnostic Use