


Product Name: Recombinant Human TENM3 Protein
Catalog #: PEH3452



Summary

Name	TEN3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P273
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DYDSEGRLTN VTFPTGVVTN LHGDMDKAIT VDISSSREE DVSITSNLSS IDSFYTMVQD QLRNSYQIGY

image



Background

Alternative Names	Teneurin-3 (Ten-3) (Protein Odd Oz/ten-m homolog 3) (Tenascin-M3) (Ten-m3) (Teneurin transmembrane protein 3)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use