


Product Name: Recombinant Human RELN Protein
Catalog #: PEH3457



Summary

Name	RELN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P78509
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GENVQFQWKQ ENLRVGEVYE ACWALDNILI INSAHRQVVL EDSLDPVDTG NWLFFPGATV KHSCQSDGNS IYFHGNEGSE FNFATTRDVD LSTEDIQEQW SEEFESQPTG WDVLGAVIGT ECGTIESGLS

image



Background

Alternative Names Reelin (EC 3.4.21.-)

Note

For Research Use Only , Not for Diagnostic Use