



Product Name: Recombinant Human DLL1 Protein
Catalog #: PEH3470



Summary

Name	DLL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O00548
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LRLQKHRPPA DPCRGETETM NNLANCQREK DISVSIIGAT QIKNTNKKAD FHGDHSADKN GFKARYPAVD YNLVQDLKGD DTAVRDAHSK RDTKCQPQGS

image



Background

Alternative Names Delta-like protein 1 (Drosophila Delta homolog 1) (Delta1) (H-Delta-1)

Note

For Research Use Only , Not for Diagnostic Use