


Product Name: Recombinant Human STX16 Protein
Catalog #: PEH3471



Summary

Name	STX16_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14662
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AIEITTQEIT QLFHRCQRAV QALPSRARAC SEQEGRLLGN VVASLAQALQ ELSTSRHAQ SGYLKRMKNR EERSQHFFDT SVPLMDDGDD NTLYHRGFTE DQLVLVEQNT LMVEEREREI RQIVQSISDL NEIFRDLGAM IVEQGTVLDR IDYNVEQSCI KTEDGLKQLH

image



Background

Alternative Names Syntaxin-16 (Syn16)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human STX16 Protein
Catalog #: PEH3471

