


Product Name: Recombinant Human BACH1 Protein
Catalog #: PEH3473



Summary

Name	BACH1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14867
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEFLENKNVQ TPQCKLRRYQ GNAKASPPLQ DSASQTYESM CLEKDAALAL PSLCPKYRKF QKAFGTDRVR TGESSVKDIH ASVQPNERSE NECLGGVPEC RDLQVMLKCD ESKLAMEPEE TKKDPASQCP TEKSEVTPFP HNSSIDPHGL YLSLLHTYD QYGDINFAGM QNTTVLTEKP LSGTDVQECT FGESQDLPLK SDLGTREDSS VASSDRSSVE

image



Background

Alternative Names Transcription regulator protein BACH1 (BTB and CNC homolog 1) (HA2303)

Note

Product Name: Recombinant Human BACH1 Protein
Catalog #: PEH3473



For Research Use Only , Not for Diagnostic Use