


Product Name: Recombinant Human PPP1R12A Protein
Catalog #: PEH3477



Summary

Name	MYPT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14974
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RQGV DIEAAR KEEERIMLRD ARQWLNSGHI NDVRHAKSGG TALHVAAAKG YTEVLKLLIQ AGYDVNIKDY DGWTPHAAA HWGKEEACRI LVDNLCDMEM VNKVGQTA FD VADEDILGYL EELQKKQNL HSEKRDKKSP LIESTANMDN NQSQKTFKNK ETLIEPEKN ASRIESLEQE KVDEEEEGKK

image



Background

Alternative Names	Protein phosphatase 1 regulatory subunit 12A (Myosin phosphatase-targeting subunit 1) (Myosin phosphatase target subunit 1) (Protein phosphatase myosin-binding subunit)
--------------------------	--

Note



Product Name: Recombinant Human PPP1R12A Protein
Catalog #: PEH3477

For Research Use Only , Not for Diagnostic Use