


Product Name: Recombinant Human HNRNPDL Protein
Catalog #: PEH3478



Summary

Name	HNRDL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O14979
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SKNQDDGKM FIGGLSWDTS KKDLEYLSR FGEVVDCTIK TDPVTGRSRG FGFVLFKDAA SVDKVLELKE HKLDGKLIDP KRAKALKGKE PPKKVFGGGL SPDTSEEQIK EYFGAFGEIE NIELPMDTKT NERRGFCFIT YTDEEPVKKL LESRYHQIGS GKCEIKVAQP KEVYRQQQQQ QKGGRGAAAG GRGGTRGRGR GQGQNWNQGF NNYYDQGYGN YNSAYGGDQN YSGYGGYDYT GYNYGNYGYG QGYADYSGQQ STYGKASRGG

image



Background

Alternative Names	Heterogeneous nuclear ribonucleoprotein D-like (hnRNP D-like) (hnRNP DL) (AU-rich element RNA-binding factor) (JKT41-binding protein) (Protein laAUF1)
--------------------------	--



Product Name: Recombinant Human HNRNPDL Protein
Catalog #: PEH3478

Note

For Research Use Only , Not for Diagnostic Use