


Product Name: Recombinant Human SPTBN2 Protein
Catalog #: PEH3479



Summary

Name	SPTN2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15020
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LGRLRDNREQ QHFLQDCHEL KLWIDEKMLT AQDVSYDEAR NLHTKWQKHQ AFMAELAANK DWLDKVDKEG RELTLEKPEL KALVSEKLRD LHRRWDELET TTQAKARSLF DANRAELFAQ SCCALESWLE SLQAQLHSDD YGKDLTSVNI LLKKQQMLEW EMVAVREKEVE AIQAQAKALA QEDQGAGEVE RTSRAVEEKF RALCQPMRER CRRLQASREQ HQFH RDVEDE

image



Background

Alternative Names	Spectrin beta chain, non-erythrocytic 2 (Beta-III spectrin) (Spinocerebellar ataxia 5 protein)
--------------------------	--

Note

Product Name: Recombinant Human SPTBN2 Protein
Catalog #: PEH3479



For Research Use Only , Not for Diagnostic Use