



Product Name: Recombinant Human FYB1 Protein
Catalog #: PEH3480



Summary

Name	FYB1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15117
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PRNIKPPFDL KSPVNEDNQD GVTHSDGAGN LDEEQDSEGE TYEDIEASKE REKKREKEEK KRLELEKKEQ KEKEKKEQEI KKKFKLTGPI QVIHLAKACC DVKGGKNELS FKQGEQIEII RITDNPEGKW LGRTARGSYG YIKTTAVEID YDSLKLKIDS LGAPSRPIED

image



Background

Alternative Names	FYN-binding protein 1 (Adhesion and degranulation promoting adaptor protein) (ADAP) (FYB-120/130) (p120/p130) (FYN-T-binding protein) (SLAP-130) (SLP-76-associated phosphoprotein)
--------------------------	---

Note

Product Name: Recombinant Human FYB1 Protein
Catalog #: PEH3480



For Research Use Only , Not for Diagnostic Use