


Product Name: Recombinant Human STX7 Protein
Catalog #: PEH3483



Summary

Name	STX7_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O15400
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PAQLAQRIS NIQKITQCSV EIQRNLNLG TPQDSPELR QLQKQKYTN QLAKETDKYI KEFGSLPTTP SEQRQRKIQK DRLVAEFTTS LTNFQKVQRQ AAEREKEFVA RVRASSRVSG SFPEDSSKER NLVSWESQTQ PQVQVQDEEI TEDDLRLIHE RESSIRQLEA DIMDINEIFK DLGMMIHEQG DVIDSIEANV ENAEVHVQQA NQQLSRAADY

image



Background

Alternative Names Syntaxin-7

Note

Product Name: Recombinant Human STX7 Protein
Catalog #: PEH3483



For Research Use Only , Not for Diagnostic Use