



Product Name: Recombinant Human RNF13 Protein
Catalog #: PEH3484



Summary

Name	RNF13_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O43567
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HRARRNRLRK DQLKKLPVHK FKKGDEYDVC AICLDEYEDG DKLRILPCSH AYHCKCVDPW LTKTKTCPV CKQKVVPSPG DSDSDTDSSQ EENEVTEHTP LLRPLASVSA QSFGLSESR SHQNMTESSD YEEDDNEDTD SSDAENEINE

image



Background

Alternative Names	E3 ubiquitin-protein ligase RNF13 (EC 2.3.2.27) (RING finger protein 13) (RING-type E3 ubiquitin transferase RNF13)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use