


Product Name: Recombinant Human NBN Protein
Catalog #: PEH3489



Summary

Name	NBN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O60934
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	REMDDVAIED EVLEQLFKDT KPELEIDVKV QKQEEDVNVR KRPRMDIETN DTFSDEAVPE SSKISQENEI GKKRELKEDS LWSAKEISNN DKLQDDSEML PKKLLLTEFR SLVIKNSTSR NPSGINDDYG QLKNFKKFKK VTYPGAGKLP HIIGGSDLIA HHARKNTELE EWLREQEMEVQ NQHAKEESLA

image



Background

Alternative Names	Nibrin (Cell cycle regulatory protein p95) (Nijmegen breakage syndrome protein 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human NBN Protein
Catalog #: PEH3489

