


Summary

Name	TADA3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75528
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QEYEFTDDPI DVPRIKND A PNRFWASVEP YCADITSEEV RTLEELLKPP EDEAEHYKIP PLGKHYSQRW AQEDLLEEKQ DGARAAAVAD KKKGLMGPLT ELDTKDVDAL LKKSEAEQHEQ PEDGCPFGAL TQRLLQALVE ENIISPMEDS PIPDMSGKES GADGASTSPR NQNKPFVPH TKSLESRIKE ELIAQGLLES EDRPAESED EVLAELRKRQ AELKALSAHN RTKKHDLRL AKEEVSRLQEL RQRVRMADNE VMDFRQKIMA ARQKKRTPK KEKDQAWKTL

image



Background

Alternative Names	Transcriptional adapter 3 (ADA3 homolog) (hADA3) (STAF54) (Transcriptional adapter 3-like) (ADA3-like protein)
--------------------------	--

Product Name: Recombinant Human TADA3 Protein
Catalog #: PEH3492



Note

For Research Use Only , Not for Diagnostic Use