


Product Name: Recombinant Human SCO1 Protein
Catalog #: PEH3494



Summary

Name	SCO1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75880
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EKLEKERQRH IGKPLLGGPF SLTTHGTGERK TDKDYLQWL LIYFGFTHCP DVCPEELEKM IQVVDEIDSI TTLPDLTPLF ISIDPERDTK EAIANYVKEF SPKLVGLTGT REEVDQVARA YRVYSPGPK DEDEDYIVDH TIIMYLIGPD GEFLDYFGQN

image



Background

Alternative Names Protein SCO1 homolog, mitochondrial

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SCO1 Protein
Catalog #: PEH3494

