


Product Name: Recombinant Human DDAH1 Protein
Catalog #: PEH3498



Summary

Name	DDAH1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O94760
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ARAERQHQLY VGVLGSKLGL QVVELPADES LPDCVFVEDV AVVCEETALI TRPGAPSRRK EVDMMKEALE KLQLNIVEMK DENATLDGGD

image



Background

Alternative Names N(G),N(G)-dimethylarginine dimethylaminohydrolase 1 (DDAH-1)
(Dimethylarginine dimethylaminohydrolase 1) (EC 3.5.3.18) (DDAH1)
(Dimethylargininase-1)

Note

For Research Use Only , Not for Diagnostic Use