


**Product Name: Recombinant Human CLSTN1 Protein**  
**Catalog #: PEH3499**



## Summary

|                              |  |
|------------------------------|--|
| <b>Name</b>                  | CSTN1_HUMAN  |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE  |
| <b>Endotoxin level</b>       | Please contact with the lab for this information   |
| <b>Tag</b>                   | C-His  |
| <b>Accession #</b>           | O94985   |
| <b>Host</b>                  | E.coli   |
| <b>Species</b>               | Human  |
| <b>Molecular Mass</b>        | 24 kDa   |
| <b>Buffer</b>                | Phosphate buffered saline  |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.   |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.                     |
| <b>Sequence</b>              | DKKVIDCLYT CKEGLDLQVL EDSGRGVQIQ AHPSQLVLTLEGEDLGELDK<br>AMQHISYLSNRQFPTPGIRR LKITSTIKCF NEATCISVPP VDGYYMVLQP<br>EETPKISLSGV HHFARAASEF ESSEGVFLFP ELRIISTITR EVEPEGDGAE<br>DPTVQESLVS EEIVHDLDTCEVTVEGEELN HEQESLEVDM ARLQQKGIEV |

## image



## Background

**Alternative Names** Calsyntenin-1 (Alcadein-alpha) (Alc-alpha) (Alzheimer-related cadherin-like protein) (Non-classical cadherin XB31alpha) [Cleaved into: Soluble Alc-alpha (SAlc-alpha); CTF1-alpha (C-terminal fragment 1-alpha)]

## Note

**Product Name: Recombinant Human CLSTN1 Protein**  
**Catalog #: PEH3499**



---

For Research Use Only , Not for Diagnostic Use