

Product Name: Recombinant Human HEXIM1 Protein
Catalog #: PEH3500



Summary

Name	HEX11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O94992
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	WGQQRQLGK KKHRRRPSKK KRHWKPYKYL TWEKKKFFDE KQSLRASRIR AEMFAKGQPV APYNTTQFLM DDHDQEEDL KTGLYSKRAA AKSDDTSDDD FMEEGGEDG GSDGMGGDGS EFLQRDFSET YERYHTESLQ NMSKQELIKE YLELEKCLSR MEDENNRLRL ESKRLGGDDA RVRELEELD RLRAENLQLL

image



Background

Alternative Names	Protein HEXIM1 (Cardiac lineage protein 1) (Estrogen down-regulated gene 1 protein) (Hexamethylene bis-acetamide-inducible protein 1) (Menage a quatre protein 1)
--------------------------	---

Note

Product Name: Recombinant Human HEXIM1 Protein
Catalog #: PEH3500



For Research Use Only , Not for Diagnostic Use