


Product Name: Recombinant Human CRTAM Protein
Catalog #: PEH3504



Summary

Name	CRTAM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95727
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KLRKAHVIWK KENEVSEHTL ESYRSRSNNE ETSSEEKNGQ SSHPMRCMNY ITKLYSEAKT KRKENVQHSK LEEKHIQVPE

image



Background

Alternative Names	Cytotoxic and regulatory T-cell molecule (Class-I MHC-restricted T-cell-associated molecule) (CD antigen CD355)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use