


**Product Name: Recombinant Human RIGI Protein**  
**Catalog #: PEH3506**



## Summary

|                              |                                                                                                                                                                                                                |
|------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | DDX58_HUMAN                                                                                                                                                                                                    |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                            |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                               |
| <b>Tag</b>                   | C-His                                                                                                                                                                                                          |
| <b>Accession #</b>           | O95786                                                                                                                                                                                                         |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                         |
| <b>Species</b>               | Human                                                                                                                                                                                                          |
| <b>Molecular Mass</b>        | 19 kDa                                                                                                                                                                                                         |
| <b>Buffer</b>                | Phosphate buffered saline                                                                                                                                                                                      |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                             |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| <b>Sequence</b>              | SWDFKKIEKL EEYRLLKRL QPEFKTRIIP TDIISDLSEC LINQECEEIL QICSTKGMMMA<br>GAEKLVECLL RSDKENWPKT LKLALEKERN KFSELWIVEK GIKDVETEDL<br>EDKMETSIDIQ IFYQEDPECQ NLENSCPS EVSDTNLYSP                                      |

## image



## Background

|                          |                                                                                                                                                                                                                                                                   |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Alternative Names</b> | Antiviral innate immune response receptor RIG-I (DEAD box protein 58) (Probable ATP-dependent RNA helicase DDX58) (EC 3.6.4.13) (RIG-I-like receptor 1) (RLR-1) (Retinoic acid-inducible gene 1 protein) (RIG-1) (Retinoic acid-inducible gene I protein) (RIG-I) |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## Note

**Product Name: Recombinant Human RIGI Protein**  
**Catalog #: PEH3506**



---

For Research Use Only , Not for Diagnostic Use