


Product Name: Recombinant Human SERPINC1 Protein
Catalog #: PEH3510



Summary

Name	ANT3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P01008
Host	E.coli
Species	Human
Molecular Mass	50 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HGSPVDICTA KPRDIPMNPM CIYRSPEKKA TEDEGSEQKI PEATNRRVWE LSKANSRFAT TFYQHLADSK NDNDNIFLSP LSISTAFAMT KLGACNDTLQ QLMEVFKFDT ISEKTSQDIH FFFAKLNCRL YRKANKSSKL VSANRLFGDK SLTFNETYQD ISELVYGAKL QPLDFKENAE QSRAAINKWV SNKTEGRITD VIPSEAINEL TVLVLVNTIY FKGLWWSKFS PENTRKELFY KADGESCSAS MMYQEGKFRY RRVAEGTQVL ELPFKGDDIT MVLILPKPEK SLAKVEKELT PEVLQEWLDE LEEMMLVVHM PRFRIEDGFS LKEQLQDMGL VDLFSPEKSK LPGIVAEGRD DLYVSDAFHK AFLEVNEEGS EAAASTAVVI AGRSLNPNRV TFKANRPFLV FIREVPLNTI IFMGRVANPC VK

image



Background

Product Name: Recombinant Human SERPINC1 Protein
Catalog #: PEH3510



Alternative Names

Antithrombin-III (ATIII) (Serpin C1)

Note

For Research Use Only , Not for Diagnostic Use