


Product Name: Recombinant Human FES Protein
Catalog #: PEH3521



Summary

Name	FES_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P07332
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LEPGELQLNE LTVESVQHTL TSVTDELAVA TEMVFRRQEM VTQLQQELRN EEENTHPRER VQLLGKRQVL QEALQGLQVA LCSQAKLQAQ QELLQTKLEH LGPGEPPPVL LLQDDRHSTS

image



Background

Alternative Names	Tyrosine-protein kinase Fes/Fps (EC 2.7.10.2) (Feline sarcoma/Fujinami avian sarcoma oncogene homolog) (Proto-oncogene c-Fes) (Proto-oncogene c-Fps) (p93c-fes)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FES Protein
Catalog #: PEH3521

