


Product Name: Recombinant Human RALB Protein
Catalog #: PEH3526



Summary

Name	RALB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P11234
Host	E.coli
Species	Human
Molecular Mass	65 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAANKSKGQS SLALHKVIMV GSGGVGKSAL TLQFMYDEFV EDYEPTKADS YRKKVVL DGE EVQIDILDTA GQEDYAAIRD NYFRSGEGFL LVFSITEHES FTATAEFREQ ILRVKAEEDK IPLL VVG NKS DLEERRQVPV EEARSKAE EW GVQYVETSAK TRANVDKVFF DLMREIRTKK MSENKDKNGK KSSKNKKSFK ERCCLL

image



Background

Alternative Names Ras-related protein Ral-B (EC 3.6.5.2)

Note

Product Name: Recombinant Human RALB Protein
Catalog #: PEH3526



For Research Use Only , Not for Diagnostic Use