


Summary

Name	NFH_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P12036
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LQALNDRFAG YIDKVRQLEA HNRSLEGEAA ALRQQQAGRS AMGELYEREV REMRGAVLRL GAARGQLRLE QEHLLEDIAH VRQLDDEAR QREEEAAAAR ALARFAQEAE AARVDLQKKA QALQEECGYL RRHHQEEVGE LLGQIQGSGA AQAQMQAETR DALKCDV TSA LREIRAQLEG HAVQSTLQSE EWFRVRLDRL SEAAKVNTDA MRSAQEEITE YRRQLQARTT ELEALKSTKD SLERQRSELE DRHQADIASY QEAIQQLDAE LRNTKWEMAA QLREYQDLLN VKMALDIEIA AYRKLLEGEE CRIGFGPIPF SLPEGLPKIP SVSTHIKVKV EEEKIKVVEKS EKETVIVVEEQ TEETQVTEEV TEEEEEKEAKE

image



Background

Alternative Names Neurofilament heavy polypeptide (NF-H) (200 kDa neurofilament protein)

Product Name: Recombinant Human NEFH Protein
Catalog #: PEH3527



(Neurofilament triplet H protein)

Note

For Research Use Only , Not for Diagnostic Use