


Product Name: Recombinant Human PTHLH Protein
Catalog #: PEH3528



Summary

Name	PTHR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P12272
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HHLIAEIHTA EIRATSEVSP NSKPSNTKN HPVRFSGSDE GRYLTQETNK VETYKEQPLK TPGKKKKGKP GKRKEQEKKK RRTRSAWLDS

image



Background

Alternative Names Parathyroid hormone-related protein (PTH-rP) (PTHrP) (Parathyroid hormone-like protein) (PLP) [Cleaved into: PTHrP[1-36]; PTHrP[38-94]; Osteostatin (PTHrP[107-139])]

Note

For Research Use Only , Not for Diagnostic Use