


Product Name: Recombinant Human CEBPB Protein
Catalog #: PEH3532



Summary

Name	CEBPB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P17676
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GAAPAPSQVK SKAKKTVDKH SDEYKIRRER NNI AVRKSRD KAKMRNLETQ HKVLELTAEN ERLQKKVEQL SRELSTLRNL FKQLPEPLLA

image



Background

Alternative Names	CCAAT/enhancer-binding protein beta (C/EBP beta) (Liver activator protein) (LAP) (Liver-enriched inhibitory protein) (LIP) (Nuclear factor NF-IL6) (Transcription factor 5) (TCF-5)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use